

## Instructions for use



PEOPLE HAVE PRIORITY



**Air motor RC-20 BC/RM**  
**Straight handpiece RC-43**  
**Contra-angle handpiece RC-58**

# Contents

---

|   |    |
|---|----|
| <b>Symbols</b> .....                      | 4  |
| in the Instructions for use .....         | 4  |
| on the medical device / packaging.....    | 5  |
| <b>1. Introduction</b> .....              | 6  |
| <b>2. Safety notes</b> .....              | 10 |
| <b>3. Product description</b> .....       | 13 |
| Air motor.....                            | 13 |
| Straight and contra-angle handpiece ..... | 14 |
| <b>4. Operation</b> .....                 | 15 |
| Assembly/Removal .....                    | 15 |
| To change rotary instruments .....        | 18 |
| Test run .....                            | 20 |
| <b>5. Hygiene and maintenance</b> .....   | 21 |
| General notes.....                        | 21 |

|   |           |
|---|-----------|
| Limitations on processing .....                     | 23        |
| Initial treatment at the point of use .....         | 24        |
| Manual cleaning .....                               | 25        |
| Manual disinfection .....                           | 27        |
| Drying .....  | 28        |
| Inspection, Maintenance and Testing .....           | 29        |
| Packaging .....                                     | 33        |
| Sterilization.....                                  | 34        |
| Storage.....  | 36        |
| <b>6. Servicing .....</b>                           | <b>37</b> |
| <b>7. W&amp;H Accessories and spare parts .....</b> | <b>38</b> |
| <b>8. Technical data .....</b>                      | <b>41</b> |
| <b>9. Disposal .....</b>                            | <b>44</b> |
| <b>Explanation of warranty terms .....</b>          | <b>45</b> |
| <b>Authorized W&amp;H service partners.....</b>     | <b>46</b> |

## Symbols

in the Instructions for use



**WARNING!**  
(risk of injury)



**ATTENTION!**  
(to prevent  
damage occurring)



General explanations,  
without risk to  
persons or objects



Do not dispose of  
with domestic waste

on the medical device / packaging



Medical Device

## Symbols

on the medical device / packaging



CE marking  
with identification number  
of the Notified Body



DataMatrix Code  
for product information  
including UDI (Unique  
Device Identification)



Data structure in  
accordance with  
Health Industry Bar Code



Catalogue number



UL Component Recognition Mark indicates compliance with  
Canadian and U.S. requirements



Serial number



Caution! According to Federal law, this medical device may only be  
sold by or on the order of a dentist, physician or any other medical  
practitioner licensed by the law of the State in which he or she  
practices and who intends to use or order the use of this medical  
device



Date of manufacture



Sterilizable up to the  
stated temperature

## 1. Introduction

---

Customer satisfaction has absolute priority in the W&H quality policy. This medical device has been developed, manufactured and subjected to final inspection according to legal regulations, quality and industry standards.

### **For your safety and the safety of your patients**

Prior to initial use please read the Instructions for use. These explain how to use your medical device and guarantee a smooth and efficient operation.



Observe the safety notes.

**Intended use (straight and contra-angle handpiece)**

The dental straight/contra-angle handpiece is intended for the following applications: Removal of decayed materials, cavities and crown preparation, removal of fillings, finishing and polishing of tooth and restoration surfaces.

**Intended use (air motor)**

The air motor is intended for the following applications: Drive for dental handpieces for dental restoration and prophylaxis. Supply of cooling air, spray air, spray liquid and light.



Misuse may damage the medical device and hence cause risks and hazards for patient, user and third parties.

### **Qualifications of the user**

We have based our development and design of the medical device on the »dentists, dental hygienists, dental employees (prophylaxis) and dental assistants« target group.

### **CE Production according to EU Directive**

0297 The medical device meets the requirements of Directive 93/42/EEC.

### **Responsibility of the manufacturer**

The manufacturer can only accept responsibility for the safety, reliability and performance of the medical device when it is used in compliance with the following directions:

- > The medical device must be used in accordance with these Instructions for use.
- > The medical device has no components that can be repaired by the user.
- > Only the components approved by the manufacturer may be replaced (Spray connection ring, Sprayclip, O-ring, coolant hose).



### **Skilled application**

The medical device is intended only for skilled application according to the intended use as well as in compliance with the valid health and safety at work regulations, the valid accident prevention regulations and in compliance with these Instructions for use.

The medical device should be prepared for use and maintained by staff who have been trained in procedures for infection control, personal safety and patient safety.

Improper use, (e.g., through poor hygiene and maintenance), non-compliance with our instructions or the use of accessories and spare parts which are not approved by W&H, invalidates all claims under warranty and any other claims.



- > Always ensure the correct operating conditions and cooling function.
- > Always ensure that sufficient and adequate cooling is delivered and ensure adequate suction.
- > In case of coolant supply failure, the medical device must be stopped immediately.
- > Use only filtered, oil-free and cooled air supplied by dental compressors to operate the medical device.
- > Check the medical device for damage and loose parts each time before using.
- > Do not operate the medical device if it is damaged.
- > Perform a test run each time before using.
- > Avoid overheating at the treatment site.
- > Do not touch the soft tissue with the head of the medical device. Risk of burning if the medical device overheats!
- > It is imperative to comply with the concentrations and exposure times specified by the manufacturer of the treatment water decontamination system, as well as its handling.



- > Before using the medical device for the first time, store it at room temperature for 24 hours.
- > The operation of the medical device is permitted only on supply units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).



### **Air motor only.**

- > Use only the supply hoses as specified by EN ISO 9168.
- > Replace faulty or leaky O-rings immediately.
- > Always follow recommendations made by the manufacturer of the transmission handpieces and the rotary instrument.

## Hygiene and maintenance prior to initial use




- > The medical device is not sterilized when delivered.
- > The packaging is non-sterilizable.



- > Clean, disinfect and lubricate the medical device.
- > Sterilize the medical device.

### 3. Product description

Air motor




#### RC-20 BC/RM


- ① Borden 2/3-hole
- ② Standard 4-hole
- ③ O-ring
- ④ Seal
- ⑤ Union nut

## Product description

## Straight and contra-angle handpiece



- ① RC-43
- ② RC-58
- ③ Chuck lever
- ④ Latch




Do not assemble or remove the medical device during operation!

- 1 Insert the air motor with BC/RM connection into the apertures of the supply hose.
- 2 Screw the union nut on.



Verify full engagement.

- 3 Check leak tightness. (not possible with BC connection, because the return air is expelled through the outer sheath)
- 4 Unscrew the union nut and remove the medical device from the supply hose.
- 5 Remove the protection cap from the external coolant tube



Do not assemble or remove the straight/contrangle handpiece during the operation!

- 1** Push the medical device onto the motor until it engages audibly. **2**



- 3** Verify full engagement

or


- 4** Remove the straight/contrangle handpiece.



## Rotary instruments



- > Use only rotary instruments which are in perfect condition. Follow the operating instructions of the manufacturer.
- > Insert the rotary instrument only when the medical device is stationary.
- > Never touch the rotary instrument while it is still rotating.
- > Do not activate the chucking system of the medical device during operation. This leads to detachment of the rotary instrument, damage to the chucking system and/or heating up of the medical device. Risk of burning!




### RC-43

- 1** Open the chuck lever.
- 2** Insert the rotary instrument up to the limit stop.
- 3** Close the chuck lever.



- 4** Verify full engagement.
- 5** Open the chuck lever and remove the rotary instrument **6**.  
Close the chuck lever.



### RC-58

- 1** Open the latch.
- 2** Insert the rotary instrument up to the limit stop and turn until it engages.
- 3** Close the latch.



- 4** Verify full engagement.
- 5** Open the latch and remove the rotary instrument **6**.  
Close the latch.

## Test run



Do not hold the medical device at eye level.

## RC-20 BC / RM

- > Start the medical device idle for 5 seconds.
- >

## RC-43, RC-58

- > Insert the rotary instrument.
- > Operate the medical device.



In the event of operating malfunctions (e.g., vibrations, unusual noise, overheating, coolant failure or leakage) **stop the medical device immediately** and contact an authorized W&H service partner.



Follow your local and national laws, directives, standards and guidelines for cleaning, disinfection and sterilization.



The information on the validated reprocessing methods serves as an example for an ISO 17664 compliant reprocessing of the medical device.



- > Wear protective clothing, safety glasses, face mask and gloves.
- > Remove the transmission instrument from the air motor.
- > Remove the air motor from the supply hose.



- > Use only oil-free, filtered compressed air with a maximum operating pressure of 3 bar for manual drying.

## Cleaning agents and disinfectants



- > Read the notes, follow the instructions and heed the warnings provided by the manufacturers of cleaning agents and/or disinfectants.
- > Use only detergents which are intended for cleaning and/or disinfecting medical devices made of metal and plastic.
- > It is imperative to comply with the concentrations and exposure times specified by the manufacturer of the disinfectant.
- > Use disinfectants which have been tested and found effective by the Verbund für Angewandte Hygiene e.V. (VAH = Association for Applied Hygiene), the Österreichischen Gesellschaft für Hygiene, Mikrobiologie und Präventivmedizin (ÖGHMP = Austrian Society for Hygiene, Microbiology and Preventive Medicine), the Food and Drug Administration (FDA) or the U.S. Environmental Protection Agency (EPA).



The user is responsible for validating its process if the specified cleaning agents and disinfectants are not available.



The product lifetime and the medical device's ability to operate correctly are mainly determined by mechanical stress during use and chemical influences due to processing.

- > Send worn or damaged medical devices and/or medical devices with material changes to an authorized W&H service partner.



### **Processing cycles**

- > Safe use is guaranteed until at least 1,000 reprocessing cycles.



Clean the medical device immediately after every treatment, to flush out any liquid (e.g., blood, saliva etc.) and to prevent settling on the internal parts.

- > Operate the medical device for at least 10 seconds at idle speed.
- > Ensure that all coolant outlets are rinsed out.



- > Wipe the entire surface of the medical device with disinfectant.
- > Remove the rotary instrument.
- > Remove the straight and contra-angle handpiece from the airmotor.
- > Remove the airmotor from the tube.




Note that the disinfectant used during pre-treatment is only for personal protection and cannot replace the disinfectant step after cleaning.






Do not place the medical device in liquid disinfectant or in an ultrasonic bath.

- > Clean the medical device under running tap water (< 35°C / < 95°F).
- > Rinse and brush off all internal and external surfaces.
- > Move moving parts back and forth several times.
- > Remove any liquid residues using compressed air.
- > When using the external coolant supply, remove the spray clip and the coolant hose.




### Clean spray nozzles

- 1 Clean coolant outlets carefully with the nozzle cleaner to remove dirt and deposits.

 The nozzle cleaner can be cleaned in an ultrasonic bath and/or in the washer-disinfector.

### Cleaning the coolant tube

- 2 Blow through the coolant tube using compressed air.

 In the case of blocked coolant outlets or coolant tubes contact an authorized W&H service partner.



W&H recommends wiping down with disinfectant.



Evidence of the medical device's basic suitability for effective manual disinfection was provided by an independent test laboratory using the »mikrozid® AF wipes« disinfectant (Schülke & Mayr GmbH, Norderstedt).




- > Ensure that the medical device is completely dry internally and externally after cleaning and disinfection.
- > Remove any liquid residues using compressed air.

### Inspection



- > Check the medical device after cleaning and disinfection for damage, visible residual soiling and surface changes.
- > Reprocess any medical devices that are still soiled.
- > Sterilize the reassembled medical device following cleaning, disinfection and lubrication.



### Lubrication Air motor



- > Lubricate the dry medical device at least once a week or after 90 minutes of use or after every internal cleaning (WD).
- > Direct the medical device downwards

### Lubrication Straight handpiece/ Contra-angle handpiece



- > Lubricate the dry medical device immediately after cleaning and/or disinfection.

### **Recommended lubrication cycles**

- > Essential after every internal cleaning
- > Before each sterilization

or

- > After 30 minutes of use or once a day

### **With W&H Service Oil F1, MD-400**

- > Follow the instructions on the oil spray can and on the packaging.

### **With W&H Assistina**

- > Follow the instructions in the Assistina Instructions for use.

## Test after lubrication



- > Direct the medical device downwards.
- > Take the medical device into operation so that excess oil can escape.





Pack the medical device and the accessories in sterilization packages that meet the following requirements:

- > The sterilization package must meet the applicable standards in respect of quality and use and must be suitable for the sterilization method.
- > The sterilization package must be large enough for the sterilization goods.
- > The filled sterilization package must not be under tension.



W&H recommends sterilization according to EN 13060, EN 285 or ANSI/AAMI ST79.



- > Read the notes, follow the instructions and heed the warnings provided by the manufacturers of steam sterilizers.
- > The program selected must be suitable for the medical device.



- > Pull off the spray clip and the coolant hose from the medical device before sterilizing.
- > Sterilize the spray clip / coolant hose and the medical device.

## Recommended sterilization procedures

Dynamic air removal cycle (prevacuum / type B; steam-flush pressure-pulse / type S)

- > 134°C (273°F) for at least 3 minutes, 132°C (270°F) for at least 4 minutes,  
121°C (250°F) for at least 15 minutes

Gravity displacement cycle (type N)

- > 134°C (273°F) for at least 10 minutes, 132°C (270°F) for at least 15 minutes  
121°C (250°F) for at least 30 minutes



Evidence of the basic suitability of the medical device for effective sterilization was provided by an independent test laboratory using the LISA 517 B17L\* steam sterilizer (W&H Sterilization S.r.l., Brusaporto [BG]), the Systec VE-150\* steam sterilizer (Systec) and the CertoClav MultiControl MC2-S09S273\*\* steam sterilizer (CertoClav GmbH, Traun).

“Dynamic-air-removal prevacuum cycle” (type B) / “Steam-flush pressure-pulse cycle” (type S): temperature 134°C [273°F] – 3 minutes\*

“Gravity-displacement cycle” (type N): temperature 121°C [250°F] – 30 minutes\*\*

\* EN 13060, EN 285, ISO 17665 / \*\* ANSI/AAMI ST55 , ANSI/AAMI ST79



- > Store sterile goods dust-free and dry.
- > The shelf life of the sterile goods depends on the storage conditions and type of packaging.

## 6. Servicing

---

### Repairs and returns

In the event of operating malfunctions immediately contact an authorized W&H service partner.



> Ensure that the medical device has been completely processed before returning it.

## 7. W&H Accessories and spare parts

---



Use only original W&H accessories and spare parts or accessories approved by W&H.

Suppliers: W&H partners


|          |                                    |
|----------|------------------------------------|
| 000301xx | W&H Assistina                      |
| 30310000 | W&H Assistina TWIN (MB-302)        |
| 10940021 | W&H Service Oil F1, MD-400 (6 pcs) |
| 02036100 | Spray cap BC/RM                    |
| 02038200 | Spray cap ISO                      |
| 11144300 | RC-43                              |
| 11245800 | RC-58                              |

## W&H Accessories and spare parts

---

|          |                      |
|----------|----------------------|
| 11042000 | RC-20 BC             |
| 11042001 | RC-20 RM             |
| 04860700 | RC-E Spray kit       |
| 02083500 | Intermediate adaptor |
| 01000700 | Seal BC              |
| 02207300 | Seal RM              |

## Spray kit



### RC-E Spray kit

- ① Sprayclip RC-43
- ② Sprayclip RC-58
- ③ Spray connection ring
- ④ 2 x coolant hose
- ⑤ 2 x O-ring
- ⑥ Nozzle cleaner



- > Perform a test run.
- > Repeat the complete hygiene and maintenance process.



## 8. Technical data

|   |  | RC-20 RM / BC              |
|---|--|----------------------------|
| Coupling  | hose-side according to standard<br>Motor/transmission instrument connection according to standard<br>Outer diameter of the motor sheath (mm) | ISO 9168<br>ISO 3964<br>18 |
| Operating pressure range  | (bar/psi)  | 2.2 – 3 bar / 32 - 43.5psi |
| Recommended operating pressure  |  | 2.5 bar / 36.3psi          |
| Speed range (rpm) at an operating pressure up to:<br>(at resultant exhaust air pressure of maximum 0.25 bar/3.6psi) |  | 25,000** +/- 10%           |
| Speed control   |  | no                         |
| Maximum torque up to  | (Ncm)  | 4**                        |
| Maximum power up to   | (W)  | 30**                       |
| max. air consumption (NI/min) at recommended operating pressure   |  | < 60                       |
| Spray water flow acc. to ISO 14457  | (ml/min)   | > 50                       |
| Water pressure  | (bar)  | 0.5 – 2*                   |
| Chip air pressure   | (bar)  | 1 – 2.5*                   |

\* Chip air pressure / water pressure must be set at the same time  
Chip air pressure must be higher than water pressure.

\*\* Power and speed are dependent on the quality of the supply hoses used and may differ from the specified values.  
rpm = min-1 (Revolutions per minute)

## Technical data

|  |                            | RC-58             | RC-43    |
|--|----------------------------|-------------------|----------|
| Coupling   | Hose-side acc. to standard |                   |          |
| Motor coupling   |                            | ISO 3964          |          |
| Speed range at 2.2 bar – 3 bar                                     | (min <sup>-1</sup> )       |                   |          |
| Maximum rated speed  | (min <sup>-1</sup> )       | 25.000 +/-10%     |          |
| Stationary torque  | (Ncm)                      |                   |          |
| Max. power   | (W)                        |                   |          |
| Air consumption at 2.2 bar – 3 bar                                 | (NI/min)                   |                   |          |
| Operating pressure   | (bar)                      |                   |          |
| Transmission ratio   |                            | 1:1               |          |
| Rotary instrument ISO 1797   | (ø mm)                     | 2,35              |          |
| Length approved by W&H   | (mm)                       | 34*               | 50*      |
| Min. chuck length  | (mm)                       | engaging          | engaging |
| Spray flow rate RC-E   | ISO 14457 (ml/min)         | > 50 (ISO 7785-2) |          |
| Recommended water pressure RC-E                                    | (bar)                      | 0,5 – 2 **        |          |
| Chip air consumption RC-E (must be higher than the water pressure) | (bar)                      |                   |          |



\* When using longer rotary instruments the user must ensure by correct selection of the operating conditions, that there is no danger to the user, patient or third parties.



### **Temperature information**

Temperature of the medical device on the operator side:

maximum 55°C (131°F)

Temperature of the medical device on the patient side:

maximum 50°C (122°F)

Temperature of the working part (rotary instrument):

maximum 41°C (105.8°F)

### **Ambient conditions**

Temperature during storage and transport:

-40°C to +70°C (-40°F to +158°F)

Humidity during storage and transport:

8% to 80% (relative), non-condensing

Temperature during operation:

+10°C to +35°C (+50°F to +95°F)

Humidity during operation:

15% to 80% (relative), non-condensing

## 9. Disposal

---



Ensure that the parts are not contaminated on disposal.



Follow your local and national laws, directives, standards and guidelines for disposal

> Medical device

> Packaging

# Explanation of warranty terms

This medical device has been manufactured with great care by highly qualified specialists. A wide variety of tests and controls guarantee faultless operation. Please note that claims under warranty can only be validated when all the directions in the Instructions for use have been followed.

**As manufacturer, W&H is liable for material or manufacturing defects within the warranty period.**

In case of complaint, please contact your nearest W&H service partner.

We accept no responsibility for damage caused by incorrect handling or by repairs carried out by third parties not authorized to do so by W&H!

Claims under warranty accompanied by proof of purchase, must be sent to the vendor or to an authorized W&H service partner. The provision of service under warranty extends neither the warranty period nor any other guarantee period.

warranty

## Authorized W&H service partners

---

Find your nearest authorized W&H service partner at <http://wh.com>  
Simply go to the menu option »Service« for full details.

Or simply scan the QR code.





## **Manufacturer**

**W&H Dentalwerk Bürmoos GmbH**

**Ignaz-Glaser-Straße 53, 5111 Bürmoos, Austria**

**t + 43 6274 6236-0,  
office@wh.com**

**f + 43 6274 6236-55  
wh.com**

**Form-Nr. 50565 AEN  
Rev. 005 / 05.05.2020  
Subject to alterations**