

SleeperOne



Comfort & Precision
*Anaesthetise with complete
peace of mind*



Dentalhitec



SleeperOne


Why choose electronic anaesthesia with SleeperOne?

Electronic injection allows perfect control of injection speeds for:

- Performing painless mucosal and palatal anaesthesia
- Optimise all infiltration, para-apical, block and intraligamentary anaesthetics.

In addition to being reliable, this system requires few consumables. It will also give a modern and differentiating image to your practice.

- Improves the effectiveness of anaesthesia
- Painless delivery thanks to electronically assisted injection
- Reduces patient anxiety, especially for children
- No increased muscular effort for dense tissue injections
- Reduces the risk of intraligamentary and palatal necrosis



	SleeperOne	Manual Syringe	Ratchet Syringe
Painless during injections	★★★★	★	★
Comfortable for children	★★★★	★	★★
Less risk of postoperative complications	★★★★	★	★
Predictability of injection comfort	★★★★	★	★★
No muscular effort	★★★★	★	★★
Fully controlled injection pressure	★★★★	★	★
Reduction of MSDs	★★★★	★	★★

Compact Kit with Wireless Foot Pedal



“ 1 ml/min is the optimal anaesthetic injection speed.

Handbook of Local Anaesthesia
Stanley F.Malamed

”

3 injection modes

Adapts to the resistance of the tissues

Aspiration function

Sterilisable cartridge container

Optimised comfort with Effitec needles thanks to the patented scalpel bevel which incises without tearing the tissue



Reassuring

The appearance reassures children and adults alike.



Precision

The pen grip places your fingers as close as possible to the needle for better stability control.



Innovation

The SleeperOne unit calculates the resistance and time, adjusting the speed and pressure accordingly.



Intuitive

The handling is easy, pleasant and does not change your protocols.

DHT Needle	30G - 16mm Ø 0.3mm	30G - 9mm Ø 0.3mm	27G - 16mm Ø 0.4mm	27G - 35mm Ø 0.4mm
Hub Colour				
Infiltration (free mucosa)	●			
Infiltration (attached gingiva)	●	●	●	
Intraseptal adult		●		
Intraseptal child		●		
PDL - Intraligamentary		●		
Palatal	●	●		
ID Block				●

New Wireless System
Now Available...

Soan



Find out more at
swallowdental.co.uk

SOAN Wireless Injection Pen

Ensure maximum comfort for you and your patient

Real pen grip with customisable injection pressure. Let the unit do the work for you with tissue resistant indicators which adapt the injection speed in real time. Use Soan's wireless, battery-free foot pedal with either single touch or constant pressure delivery options.

My DHT offers 360° support

- Easy access to all your services (training, device, videos, etc).
- A complete guide to getting started with SleeperOne.
- Clinical videos by sector, for adults and children.
- Clinical question?
Our experienced clinicians are ready to answer.



☎ 01535 656 312 ✈ [SWALLOWDENTAL.CO.UK](http://swallowdental.co.uk)

✉ SALES@SWALLOWDENTAL.CO.UK

SWALLOW